



## Protocol for Filter Aided Sample Preparation (FASP)

The use of detergents to lyse cells can be problematic for down-stream mass spectrometry. Routinely used detergents such as Triton-X and NP-40 can cause ion suppression in the mass spectrometer while SDS can additionally impair enzyme digestion. Filter aided sample preparation or FASP is a method for generating tryptic peptides (for LC-MS/MS) from detergent lysed cells and tissues, that incorporates a buffer exchange step for the removal of the detergent. (Matthias Mann, <http://www.biochem.mpg.de/226356/FASP>).

### *MATERIALS*

#### NOTES

- The following items of electrical equipment are used in this procedure. You should ensure that you have been correctly instructed in their correct use prior to carrying out this procedure.
  - Bench-top microfuge
  - Heating Block
  - Sonicator
  - Speedivac
  - Thermomixer
  - Water Bath
- This protocol involves an overnight step
- Use Eppendorf® LoBind microfuge tubes Protein
- Filter units - Sartorius, Vivacon 500, Product number: VN0H22

#### Buffers

Note: Urea and Iodoacetamide buffers must be *freshly prepared*.  
Use HPLC grade water.

Lysis Buffer 0.1% SDS; 50mM-Tris pH 7.5

Urea Buffer 8M-Urea; 0.1M-Tris pH 8.5 (4.8 g Urea per 10 ml, 1 ml 1M-Tris pH 8.5) *fresh*

IAA Buffer 0.05M-Iodoacetamide (5 mg in 500  $\mu$ L 8M-Urea/Tris buffer) *fresh*

Ambic 0.05M-Ammonium bicarbonate (0.4 g in 100 ml water)



## **METHOD**

### **Sample lysis**

Methods of sample lysis other than the following are available. The biological sample (for example a single 50  $\mu\text{L}$  cell pellet from a 15 cm dish) is first lysed in approximately 150  $\mu\text{L}$  of a buffer containing SDS (1% SDS; 50 mM-Tris pH 7.5).

1. Incubate the sample at 95°C for 5 minutes.
2. Sonicate the sample on ice x 3 briefly to reduce viscosity and shear DNA
3. Clarify the lysate by centrifugation at 16,000 x g for 5 minutes
4. Estimate the protein content. In the following example 200  $\mu\text{g}$  of protein in 100  $\mu\text{L}$  was processed.

### **Reduction of sample**

1. Add 1M-DTT to give a final DTT concentration of 0.1M in each sample
2. Incubate on Thermomixer at 90°C for 10 minutes
3. Spin briefly

### **Buffer exchange (Urea displaces SDS)**

1. Take up to 100  $\mu\text{L}$  of sample and bring volume up to 0.5 mL with Urea buffer (8M-Urea; 0.1M-Tris pH 8.5)
2. Load onto a 10 kDa filter unit
3. Spin at 14,000 x g for 15 mins
4. Discard the flow through
5. Add 200  $\mu\text{L}$  of Urea buffer to filter
6. Spin at 14,000 x g for 15 mins
7. Discard the flow through

### **Alkylation of sample**

1. Add 100  $\mu\text{L}$  IAA buffer (0.05M Iodoacetamide in Urea Buffer) to filter and mix in Thermomixer at 600 rpm for 1 minute
2. Incubate for 30 minutes at room temperature in the dark
3. Spin at 14,000 x g for 15 mins
4. Discard flow through

### **Removal of IAA**

1. Add 200  $\mu\text{L}$  Urea buffer to the filter
2. Spin at 14,000 x g for 15 mins
3. Repeat the last three steps

### **Removal of Urea before trypsin digestion**

1. Add 200  $\mu\text{L}$  Ambic (0.05 M-Ammonium bicarbonate) to the filter
2. Spin at 14,000 x g for 15 mins



3. Discard flow through
4. Repeat the last two steps

### Trypsin digestion (overnight step)

1. For every 100  $\mu\text{g}$  of protein add 2  $\mu\text{g}$  Trypsin in a total volume of 40  $\mu\text{L}$  Ambic
2. Incubate filter units in water bath at 37°C overnight (units are too high to fit in thermomixer)
3. Transfer filter units to new collection tubes
4. Centrifuge the filters at 14,000 x g for 20 mins
5. RETAIN FILTRATE – THIS IS THE TRYPTIC PEPTIDE MIXTURE
6. Add 40  $\mu\text{L}$  Ambic to filter
7. Spin at 14,000 x g for 15 mins
8. Transfer the collected peptide digest to clean Lo-bind tube
9. Determine concentration with Nanodrop at 280 nm (using Ambic as blank)
10. Dry the digest under vacuum
11. Resuspend the digest in 0.5% TFA to 2  $\mu\text{g}/\mu\text{L}$

*The re-suspended tryptic digest should be de-salted using Zip tips as per protocol*